YCHROMSTRIP

AZF microdeletions detection of Y chromosome

Procedure validated for automated and manual use

 **opegen**
by operon

www.opegen.com - sales@opegen.com - Tel.: +34 976 503 597 - Spain

AZF microdeletions detection of Y chromosome

The genetic basis of this disorder was discovered in 1976 by Tiepolo and Zuffardi. They found microscopic deletions in the distal portion of the Y chromosome in four men with azoospermia. On this basis, they postulated the existence of one or more critical genes for spermatogenesis in this chromosomal region which, later, was subsequently validated and referred to as Azoospermia Factor (AZF).

These microdeletions were identified in three regions of Yq11: AZFa, AZFb, and AZFc. It has shown that deletion of the AZFc region is the most commonly defined molecular cause of spermatogenic failure (from 7 to 20% according to some authors).

Intended use

YchromStrip is a test for the detection of deletions in the azoospermia factor (AZF) located on chromosome Y. These deletions are the cause of male infertility due to azoospermia (no sperm) or oligospermia (low sperm count).

The test detects the presence or absence of six STSs (sequence tagged sites) located in three different regions of Azoospermia Factor: AZFa (sY84 and sY86), AZFb (sY127 and sY134) and AZFc (sY254 and sY255). In principle, the analysis of a single non-polymorphic STS for each region would be sufficient to determine if there is a deletion; however, analysing two from each region, it is increased the accuracy of the test.

In the test also include probes to detect the presence of three genes:

- ZFY (Zinc Finger Y region) and ZFX (Zinc Finger X region), two highly homologous genes located on Y and X chromosome, respectively. They are included to control the presence of DNA and to check that there has been a proper amplification of the sample.
- SRY (Sex reverse Y), located on the short arm of Y chromosome and testicular main determinant.

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Step 3 Paso 3

Hybridation protocol
Protocolo de hibridación

* OpeGen kits includes all the necessary reagents even Taq polymerase

RESULTS Resultados

- 1) Simoni M. et al. International Journal of Andrology (2004). Vol 27: 240-249. "EAA/EMQN best practice guidelines for molecular diagnosis of Y-chromosomal microdeletions. State or the art 2004".
- 2) Vogt P. Human Reproduction Update (2005). Vol 11 (4): 319-336. "AZF deletions and Y chromosomal haplogroups: history and update based on sequence".
- 3) Kühnert B. et al. Human Reproduction (2004). Vol 19 (4): 886-888. "Case report: natural transmission of an AZFc Y-chromosomal micro-deletion from father to his sons".