MOLECULAR DIAGNOSTICS

HYBRIDIZATION ON STRIP / OLIGOIMMUNOCHROMATOGRAPHY

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
GENERAL PROCEDURE

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
* OpeGen kits include all the necessary reagents even Taq polymerase

HemochromaStrip

HemochromaStrip is a test designed for the detection of three point mutations/genetic variants of HFE gene associated with hemochromatosis: C828Y, H63D and S65C.

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
LactoStrip

LactoStrip is a test designed for the detection of two MCM6 gene polymorphisms associated with the persistence of lactase activity in adults: C/T13910 and G/A22018

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
ThromboStrip


ThromboStrip is a test designed for the detection of three point mutations/genetic variants associated with the risk of suffering deep venous thrombosis:

- G1691A mutation in Factor V;
- G20210A mutation in prothrombin; and
- C677T mutation in methylenetetrahydrofolate reductase.

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AZF microdeletions detection of Y chromosome

The genetic basis of this disorder was discovered in 1976 by Tiepolo and Zuffardi. They found microscopic deletions in the distal portion of the Y chromosome in four men with azoospermia. On this basis, they postulated the existence of one or more critical genes for spermatogenesis in this chromosomal region which, later, was subsequently validated and referred to as Azoospermia Factor (AZF).

These microdeletions were identified in three regions of Yq11: AZFa, AZFb, and AZFc. It has shown that deletion of the AZFc region is the most commonly defined molecular cause of spermatogenic failure (from 7 to 20% according to some authors).

Intended use

YchromStrip is a test for the detection of deletions in the azoospermia factor (AZF) located on chromosome Y. These deletions are the cause of male infertility due to azoospermia (no sperm) or oligospermia (low sperm count).

The test detects the presence or absence of six STSs (sequence tagged sites) located in three different regions of Azoospermia Factor: AZFa (sY84 and sY86), AZFb (sY127 and sY134) and AZFc (sY254 and sY255). In principle, the analysis of a single non-polymorphic STS for each region would be sufficient to determine if there is a deletion; however, analysing two from each region, it is increased the accuracy of the test.

In the test also include probes to detect the presence of three genes:

- ZFY (Zinc Finger Y region) and ZFX (Zinc Finger X region), two highly homologous genes located on Y and X chromosome, respectively. They are included to control the presence of DNA and to check that there has been a proper amplification of the sample.
- SRY (Sex reverse Y), located on the short arm of Y chromosome and testicular main determinant.

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


Detection of HLA haplotypes associated to celiac disease

Celiac disease is a chronic inflammatory condition of the intestine that is triggered by the consumption of gluten or associated proteins found in wheat, barley and rye. It is one of the most common conditions affecting the Caucasian population, with a prevalence of between 1:100 and 1:500 in Europe and North America.

Susceptibility to gluten sensitivity is to some extent genetically predetermined. In most of the human populations studied, 90-95% of patients were carriers of the heterodimer HLA-DQ2, encoded by alleles DQA1*05 and DQB1*02 in cis position (more common in Central and Northern Europe) or in trans position (more common in Mediterranean countries). The remaining patients (5-10%) often have a second heterodimer, HLA-DQ8 (mainly in indigenous South American patients), encoded by alleles DQA1*03 and DQB1*0302. Patients who are non-carriers of DQ2 or DQ8 may show at least one separate DQ2 allele; very few cases have been described where both alleles are absent.

The following diagram lists the HLA haplotypes that constitute antigens HLA-DQ2 and HLA-DQ8 (the alleles DQA1 and DQB1 that encode the antigen's α and β chains, and their associated DRB1 are also indicated):

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Intended use

CeliacStrip is an assay that allows the identification of the main HLA haplotypes associated with celiac disease.

The test detects the presence or absence of haplotypes that encode HLA-DQ2 (cis and trans) and HLA-DQ8.

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
1. Petronzelly F. et al. Ann. Hum. Genet. (1997); 61:307-317. "Genetic contribution of the HLA region to the familial clustering of coeliac disease".
2. Louka A.S. et al. Tissue Antigens (2003); 61: 105-117. "HLA in coeliac disease: Unravelling the complex genetics of a complex disorder".
3. Kipatány A. et al. Aliment. Pharmacol. Ther. (2006); 24 (9): 1395-1402. "Diagnostic significance of HA-DQ typing in patients with previous coeliac disease diagnosis based on histology alone".
4. Ferrer A et al. Biotecnología aplicada (2005); 22: 91-101. "Overview on HLA and DNA typing methods".

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Results interpretation


The presence of some of the haplotypes associated with celiac disease will be determined by the presence of the group of bands that identify the alleles that form each one of the haplotypes. Therefore:

- A sample will be positive for DQ2 cis when at least the bands corresponding to haplotypes DQA1*05 - DQB1*02 - DRB1*03 are present, in addition to the β -globin (BGB) band.
- A sample will be positive for DQ2 trans when at least the bands corresponding to haplotypes DQA1*05 - DQB1*0301 - DRB1*11 o DRB1*12 and DQA1*02 - DQB1*02 - DRB1*07 are present, in addition to the β -globin (BGB) band.
- A sample will be positive for DQ8 when at least the bands corresponding to haplotype DQA1*03 - DQB1*0302 - DRB1*04 appear.
- A sample that does not carry one of the haplotypes associated with celiac disease could yield:
 - the β -globin (BGB) band only, or
 - some of the bands corresponding to the alleles that constitute the haplotypes detected with the kit, but in the absence of others (incomplete haplotypes).

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Detection of HLA B27 alleles

The existence of a genetic predisposition to various diseases in individuals carrying the HLA-B*27 histocompatibility antigen has been demonstrated: ankylosing spondylitis, acute anterior uveitis, reactive arthritis, inflammatory bowel disease.

Ankylosing spondylitis (or ankylosing spondylarthritis, AS) is a chronic, disabling rheumatic disease, more common especially in white people (0.5-1% of the population). It usually appears in males between 20 and 30 years of age.

Acute anterior uveitis is the most common form of endogenous uveitis. Uveitis is an inflammation of the uveal tract, a pigmented vascular layer that lies between the sclera and the retina. The main symptoms are photophobia, pain, erythema, decreased vision and tearing.

Reactive arthritis is an inflammation of one or more joints that appears as a reaction to an infection and which usually causes inflammation of the urinary and genital tract, the joints and eyes.

Inflammatory bowel disease refers to a group of disorders that predominantly affect the intestine and which are characterized by chronic inflammation, which tends not to subside. It brings together various diseases, but especially Crohn's disease and ulcerative colitis. The existence of a genetic predisposition to various diseases in individuals carrying the HLA-B*27 histocompatibility antigen has been demonstrated: ankylosing spondylitis, acute anterior uveitis, reactive arthritis, inflammatory bowel disease.

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Intended use


HLA B27 Strip test allows the determination of the HLA B*27 allele. This group of alleles has been associated with a wide range of pathologies such as ankylosing spondylitis (88-95% of patients), acute anterior uveitis (50-60% of patients), reactive arthritis (60-85% of patients), and inflammatory bowel disease (50-60% of patients).

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
Detection of HLA B5701 alleles

Abacavir is used in the treatment of HIV-1 infection since 1999. The WHO recommends abacavir as a second-line treatment, owing to the risk of hypersensitivity associated with its use, which affects 5-8% of Caucasian patients.

Several studies have demonstrated the existence of a strong predictive correlation, in Caucasian populations and groups of Hispanic ethnicity, between hypersensitivity to abacavir and the presence of the HLA-B*5701 allele. This correlation is strong enough to enable us to predict the risk of hypersensitivity to abacavir and classify individuals as low (<1%) or high (>70%) risk, based on the absence or presence of the HLA-B*5701 allele.

Several polymorphisms has been identified. These ones explain nearly 15% of the variability in HIV-1 viral load in individuals during the asymptomatic phase of the disease were identified. One of them, located in the HCP5 gene (HLA complex P5), shows a perfect linkage disequilibrium with HLA-B*5701 ($r^2 = 1$).

The useful role of this SNP in predicting hypersensitivity to abacavir has been demonstrated in European populations (100% detection sensitivity for HLA-B*5701, specificity of 99.4%, positive predictive value of 94.2% and negative predictive value of 100%).

Intended use

HLA B5701 Strip is a test designed for the detection of HLA B5701 alleles associated with Abacavir hypersensitivity reaction.

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RESULTS


Resultados

A B*5701 positive sample must always be positive for the bands S1 and S2 as well as for the bands associated with the two control genes (Control A and Control B).

All results in which the bands corresponding to the two amplification controls fail to appear must be regarded as invalid.

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
- 1) Saag M et al. "High sensitivity of human leukocyte antigen-B*5701 as a marker for immunologically confirmed abacavir hypersensitivity in white and black patients". *Clinical Infectious Diseases* (2008); 46: 1111-1118.
- 2) Mallal M.B. et al. "HLA-B*5701 screening for hypersensitivity to Abacavir". *The New England Journal of Medicine* (2008); 358 (6): 568-579
- 3) Colombo S. et al. "The HCP5 single-nucleotide polymorphism: a simple screening tool for prediction of hypersensitivity to reaction to abacavir". *JID* (2008); 198: 864-867.

Introduction

Human *Papillomavirus* (HPV) is responsible for a variety of skin and mucosa disorders. To date, over 100 different HPV genotypes have been described, of which more than 40 infect the anogenital tract. Of these, approximately 1/3 is associated with cervical cancer and anal neoplasm. The anogenital HPVs are generally divided into two categories: those with a low oncogenic potential risk (Low-Risk group) and those with a medium-high oncogenic potential risk (High-Risk group). The high-risk HPVs are generally associated with high-grade precancerous lesions and invasive cancer, while low-risk HPVs are frequently found in asymptomatic or benign conditions such as genital warts.

Intended use

High PapillomaStrip test allows qualitative detection of the DNA from 19 human *Papillomavirus* subtypes of medium-high risk: 16, 18, 26, 31, 33, 35, 39, 45, 51, 52, 53, 56, 58, 59, 66, 68, 69, 73 and 82.

On the other hand, Low PapillomaStrip test allows the genotyping of DNA samples from 18 human *Papillomavirus* subtypes of low risk: 6, 11, 40, 42, 43, 44, 54, 61, 62, 67, 70, 71, 72, 74, 81, 83, 84 y 91.

Specific primers for E6-E7 regions for each 37 HPVs detected with the kit are used to perform the multiplex PCR. Since the expression of the E6-E7 genes is directly associated with carcinogenesis, detection of these targets may provide significant biological and pathological information. Some studies showed that in cervical cancer samples that were positive for the E6 region, the L1 gene could not be detected in approximately 30% of cases, due to a gene deletion process that occurs when integrating DNA of high-risk HPV in epithelial cells.

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

Step 2 Paso 2

PCR*

* OpeGen kits include all the necessary reagents even Taq polymerase

- ✓ Previous hybridization to the strip detection: it increases specificity and reproducibility of the process.
- ✓ All the reagents Ready To Use.
- ✓ Minimum necessary equipment.
- ✓ Kits include all the PCR reagents, even the Taq polymerase.

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



Productos

CT OligoGen

CT OligoGen is a test designed for the detection of *Chlamydia trachomatis* by simultaneous amplification of two independent targets of its DNA: one in the cryptic plasmid and another in its own genome and subsequent identification by immunochromatography. This test can detect infections caused by wild-type CT, the Swedish variant or those that do not carry the cryptic plasmid.

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
HSV1/2 OligoGen

HSV1/2 OligoGen is a test designed for the detection and differentiation of *Herpes simplex virus* type 1 and 2, by a common amplification and subsequent specific identification in two different immunochromatographic strips.

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NG OligoGen

NG OligoGen is a test designed for the detection of *Neisseria gonorrhoeae* by simultaneous amplification of two independent targets of its DNA and subsequent identification in immunochromatographic strip.

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OPERON S.A.